

#### **Tables for Food Texture Analyzer**

# **Optional tables 3**

## (Container holder table, Stick holder)

Tables to hold samples are available.

## **Container holder table (FR-BG-60)**



		i calu
It holds a	container not to slip	down.

#### Example of use

In the case where semi-solid sample in a container is measured.

#### [Measurement examples]

Viscosity of soup	
MADA TON PROMISE TONAL STATE OF THE PROMISE TONA	



1. Disk probe

1. Cylindrical probe

#### [Specifications]


Model	FR-BG-60	
Maximum diameter	m14~ m20	
of a container ØA	φ14~φ80	
Maximum height	15~100	
of a container		
Weight	295g	

<sup>\*</sup>Unit: mm

<sup>\*</sup> We also offer customized sizes of all attachments.




#### [Dimensions]



<sup>\*</sup>Unit: mm

## Stick holder (FR-ST-25)



#### [Measurement example]



1. Wedge probe (FR-KS series)

#### [Specifications]

Model	FR-ST-25	
Diameter of a sample Ø A	φ14~φ25	
Material	Aluminum	
Weight	253g	


<sup>\*</sup>Unit: mm

<sup>\*</sup>We also offer customized sizes of all attachments.

<sup>\*</sup>We also offer customized sizes of all attachments.



#### [Dimensions]



<sup>\*</sup>Unit: mm

### [Cautions]

- Information in this document is subject to change without prior notice.
- This document is product descriptions and handling precautions, and do not guarantee various characteristics or safety.
- This product is designed for force measurement purpose only.
- Do not copy and use this content without authorization.
- Some samples may not be suitable to measure with this product.

#### **IMADA CO., LTD**

99 Jinnoshinden-cho aza Kanowari Toyohashi

Aichi Japan 441-8077 Tel: +81-(0)532-33-3288 Fax: +81-(0)532-33-3866 E-mail: info@forcegauge.net

Website: http://www.forcegauge.net/en/



Visit our website for more information on a wide range of product specifications, measurement applications and videos.

<sup>\*</sup> We also offer customized sizes of all attachments.