

## Guide Ruler for Adhesive Tape Peel Test GRP-50

- Designed to help correctly position the Adhesive Tapes for the Peel Tests with the Peel Test fixtures
- With the guiding, the tapes are fixed parallel to the peeling direction for accurate measurements
- The Sample positionings are easier for repeated tests to help with the testing efficiency



# The Guide Ruler is to be used with the Peel Test Fixtures to help correctly position the Adhesive Tapes for the Peel Tests. Adjust the ruler according to the sample width Adjust the ruler according to the sample width Ruler to achieve accurate parallel positioning (for Peel Tests).




Specifications		
Model	GRP-50	
Compatible products	P90-200N, P90-200N-EZ, P180-200N, SSP-1608	
Compatible sample width	Max. 50mm	
Dimensions	See [Dimensions]	
Weight	Approx. 300g	

Related Products		
90 degree Peel Test Fixture  P90-200N	180 degree Peel Test Fixture P180-200N	
Sliding table for 90 degree peel test for adhesive tapes. Different sheet types can be pasted to the board for testing.	Testing fixture for 180 degree peel test for different types of adhesive tapes.	
Rubber Roller for Peel Test  APR-97	Film Cutter FSC-1525	
Roller rolled on the adhesive tape samples to apply even pressure for the peel tests.	Cutter for making testing paper or film in strips of specific widths.	



## [Dimensions]



Unit: mm

## [Cautions]

- Information in this document is subject to change without prior notice.
- This document is product descriptions and handling precautions, and do not guarantee various characteristics or safety.
- This product is designed for force measurement purpose only.
- Do not copy and use this content without authorization.

### **IMADA CO., LTD**

99 Jinnoshinden-cho aza Kanowari Toyohashi Japan 441-8077

Tel: +81-(0)532-33-3288 Fax: +81-(0)532-33-3866 E-mail: <u>info@forcegauge.net</u>

Website: https://www.forcegauge.net/en/



Visit our website for more information on wide product specifications, measurement applications and videos.